



## General features

- particle-filtering folding mask
- individually packed
- colour: white

## Safety features

- certified FFP2 face mask to EN 149:2001+A1:2009
- D label (dolomite dust test passed): maximum breathability even with really high dust levels
- cannot be reused

## Comfort features

- flexible and adjustable nose clip
- soft sealing lip guarantees a secure, comfortable sealed position of the mask
- soft material edges for high wearing comfort
- adjustable head band offering easy handling and comfortable fit





## FFP2 DATA SHEET AND IMAGES

**Particulate Respirator FFP2 Model: MASK-FFP2**



## WARNING

THIS RESPIRATOR HELPS PROTECT AGAINST CERTAIN PARTICLES. Misuse may result in **sickness or death**. For correct use, consult supervisor and these *User Instructions*.

### **Importer:**

Qvis Lighting and Security LTD

**Address:** 36 New Lane, Havant, Hampshire, PO9 2JL, United Kingdom

**Contact:** +44 (0)23 92488300

### **IMPORTANT**

Before use, wearer must read and understand these *User Instructions*. Please keep these instructions for future reference.

### **Use For**

Particles such as those from grinding, sanding, sweeping, sawing, bagging, or processing minerals, coal, iron ore, flour, metal, wood, pollen, and certain other substances. Liquid or non-oil-based particles from sprays that do not also emit oil aerosols or vapours. Follow all applicable local regulations.

### **Do Not Use For**

Do not use for gases and vapours, oil aerosols, asbestos, or sandblasting; particulate concentrations that exceed either 10 times the occupational exposure limit or applicable government regulations, whichever is lower. This respirator does not supply oxygen.

### **Biological Particles**

This respirator can help reduce inhalation exposures to certain airborne biological particles (e.g. mould, Bacillus anthracis, Mycobacterium tuberculosis, etc.) but cannot eliminate the risk of contracting infection, illness or disease.



## FFP2 DATA SHEET AND IMAGES

### **Instructions on how to use**

1. Failure to follow all instructions and limitations on the use of this respirator and/or failure to wear this respirator during all times of exposure can reduce respirator effectiveness and may result in sickness or death.
2. Follow all applicable local regulations.
3. The particles which can be dangerous to your health include those so small that you cannot see them.
4. Leave the contaminated area immediately and contact supervisor if dizziness, irritation, or other distress occurs.
5. Store the respirator away from contaminated areas when not in use.
6. Inspect respirator before each use to ensure that it is in good operating condition. Examine all the respirator parts for signs of damage including the two headbands, attachment points, nose foam, and nose clip. The respirator should be disposed of immediately upon observation of damaged or missing parts. Filtering facepieces are to be inspected prior to each use to assure there are no holes in the breathing zone other than the punctures around staples and no damage has occurred. Enlarged holes resulting from ripped or torn filter material around staple punctures are considered damage. Immediately replace respirator if damaged.
7. Conduct a user seal check before each use as specified in the Fitting Instructions section. If you cannot achieve a proper seal, do not use the respirator.
8. Dispose of used product in accordance with applicable regulations.

### **Use Limitations**

1. This respirator does not supply oxygen. Do not use in atmospheres containing less than 19.5% oxygen.
2. Do not use when concentrations of contaminants are immediately dangerous to life and health, are unknown or when concentrations exceed 10 times the permissible exposure limit (PEL) or according to specific applicable government regulations, whichever is lower.
3. Do not alter, wash, abuse or misuse this respirator.
4. Do not use with beards or other facial hair or other conditions that prevent a good seal between the face and the sealing surface of the respirator.
5. Respirators can help protect your lungs against certain airborne contaminants. They will not prevent entry through other routes such as the skin, which would require additional personal protective equipment (PPE).
6. This respirator is designed for occupational/professional use by adults who are properly trained in its use and limitations. This respirator is not designed to be used by children.
7. Individuals with a compromised respiratory system, such as asthma or emphysema, should consult a physician and must complete a medical evaluation prior to use.



## FFP2 DATA SHEET AND IMAGES

8. When stored in accordance with temperature and humidity conditions specified below, the product may be used until the “use by” date specified on the packaging


### **Storage Conditions and Shelf Life**

Before use, store respirators in the original packaging away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture and damaging chemicals. When stored in accordance with temperature and humidity conditions specified below, the product may be used until the “use by” date specified on packaging. Always inspect product and conduct a user seal check before use as specified in the User Instructions. If you cannot achieve a proper seal, do not use the respirator.



#### **End of Shelf Life**

Use respirators before the “use by” date specified on packaging



#### **Storage Temperature Range**

-20°C (-4°F) to +30°C (+86°F).



#### **Storage Maximum Relative Humidity**

<80% RH

### **Fitting Instructions (Must be followed each time respirator is worn).**

1. Pre-stretch both straps before placing respirator on the face (Fig. 1).
2. Use both hands to pull ear-bands, so that the respirator is in place against your chin, covering both your mouth and nose. Ensure the nosepiece is up. (Fig.2).
3. Position the ear-bands around your ears and adjust to ensure firm hold. (Fig. 3).
4. Place your fingertips from both hands at the top of the metal nosepiece. Using two hands, mould the nose area to the shape of your nose by pushing inward while moving your fingertips down both sides of the nosepiece (Fig. 4).



Pinching the nosepiece using one hand may result in improper fit and less effective respirator performance. Use two hands.

5. Perform a User Seal Check prior to each wearing. To check the respirator-to-face seal, place both hands completely over the respirator and exhale sharply. Be careful not to disturb the position of the respirator. If air leaks around nose, readjust the nosepiece as described in step 4. If air leaks at the respirator edges, work the straps back along the sides of your head (Fig. 5). If you CANNOT achieve a proper seal, DO NOT enter the contaminated area. See your supervisor.



## FFP2 DATA SHEET AND IMAGES



Removal Instructions See step 3 of Fitting Instructions and cup respirator in hand to maintain position on face. Pull one ear-strap away from ear. Still holding respirator in position, pull the second strap off the other ear.

This respirator contains no components made from natural rubber latex.

### **Cautions and Limitations:**

- 1) NOT FOR USE IN ATMOSPHERES CONTAINING LESS THAN 19.5 PERCENT OXYGEN.
- 2) NOT FOR USE IN ATMOSPHERES IMMEDIATELY DANGEROUS TO LIFE OR HEALTH.
- 3) DO NOT EXCEED MAXIMUM USE CONCENTRATIONS ESTABLISHED BY REGULATORY STANDARDS.
- 4) FAILURE TO PROPERLY USE AND MAINTAIN THIS PRODUCT COULD RESULT IN INJURY OR DEATH.
- 5) ALL APPROVED RESPIRATORS SHALL SELECTED, FITTED, USED AND MAINTAINED IN ACCORDANCE WITH LOCAL REGULATIONS.
- 6) NEVER SUBSTITUTE, MODIFY, ADD, OR OMIT PARTS.
- 7) REFER TO THESE USER INSTRUCTIONS FOR USE OF RESPIRATOR.